





# 3mL Syringe Instructions

1. Firstly, 1mL (or milliliter) is exactly equal to 1cc (cubic centimeter). The two terms can be used interchangeably.
2. The total volume of the syringe is 3mL. Therefore, it is measured in six 0.5mL increments. See #1 and #2 below.
3. Each small line in between a large line marked with a number is equal to 0.1mL. See #3 and #4 below. For example, you could have 0.2mL, 1.2mL, or 2.7mL.
4. Pull the plunger back so the tip of the black end is level with the line you are measuring to.



#1  
3mL  
syringe



#2  
0.5 mL



#3  
0.2 mL



#4  
1.4 mL