





# PZI Insulin Syringe Instructions

1. Firstly, 1u (unit) is **NOT** the same as 1mL (milliliter). A unit of insulin is *much* smaller than a milliliter (giving 1mL of insulin would probably result in death in a cat).
2. There are different types of insulin syringes depending on what type of insulin your cat is receiving. Make sure the syringe you are using is made to measure the specific type of insulin you are giving (Humulin insulin vs. PZI insulin).
3. This is a 20 unit syringe – see #1 below. Therefore, it is measured into four 5 unit increments.
4. Each small line is equal to 1 unit of insulin. See #2 and #3 below. For example, you can have 1 unit or 4 units. The first line closest to the needle equals 0 – it's actually the 2<sup>nd</sup> line on the syringe that equals 1 unit.
5. Pull the plunger back so the tip of the black end is level with the line you are measuring to.



#1  
20 unit syringe



#2  
1 unit



#3  
3 units